


GPR37 Rabbit mAb

ABZ2975

Product Overview

Name	GPR37 Rabbit mAb
Accession(Primary)	O15354
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O15354)

Gene Name GPR37**About Gene Symbol**

Western blot analysis of GPCR GPR37 in rat Brain lysates using GPR37 antibody.

