

# Growth Hormone Receptor Rabbit mAb


ABZ2977

## Product Overview

|                             |                                                          |
|-----------------------------|----------------------------------------------------------|
| <b>Name</b>                 | Growth Hormone Receptor Rabbit mAb                       |
| <b>Accession(Primary)</b>   | P10912                                                   |
| <b>Host Species</b>         | Rabbit                                                   |
| <b>Clonality</b>            | Monoclonal                                               |
| <b>Isotype</b>              | IgG                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                        |
| <b>Calculated MW</b>        | Calculated MW: 72 kDa; Observed MW: 120 kDa              |
| <b>Recommended Dilution</b> | WB-1:1000-1:15000; ICC/IF-1:100-1:200;<br>FC-1:200-1:500 |

## Target information (P10912)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | GHR |
| <b>About Gene Symbol</b> |     |



Western blot analysis of Growth hormone receptor in rat Brain lysates using Growth Hormone Receptor antibody.