

GRP78 BiP Rabbit mAb


ABZ2978

Product Overview


Name	GRP78 BiP Rabbit mAb
Accession(Primary)	P11021
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (P11021)

Gene Name	HSPA5
About Gene Symbol	HSPA5 (Heat Shock Protein Family A (Hsp70) Member 5) is a molecular chaperone involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mecha



Western blot analysis of GRP78 BiP in C6 lysates using GRP78 BiP antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GRP78 BiP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.