

GST3 / GST pi Rabbit mAb


ABZ2979

Product Overview


Name	GST3 / GST pi Rabbit mAb
Accession(Primary)	P09211
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100; IF-1:20-1:50; FC-1:20-1:50

Target information (P09211)


Gene Name	GSTP1
About Gene Symbol	GSTP1 (Glutathione S-Transferase Pi 1) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in hearing loss and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of GST3 (green) in A549 using GST3 antibody, and DAPI (blue).



Western blot analysis of GST3/GST pi in K562, C6, 3T3, HeLa lysates using GST3/GST pi antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using GST3/GST pi antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.