


HDAC11 Rabbit mAb

ABZ2982

Product Overview

Name	HDAC11 Rabbit mAb
Accession(Primary)	Q96DB2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q96DB2)

Gene Name HDAC11**About Gene Symbol**

Western blot analysis of HDAC11/HD11 in T47D, MCF-7 lysates using HDAC11 antibody.

