

HDAC3 Rabbit mAb


ABZ2983

Product Overview

Name	HDAC3 Rabbit mAb
Accession(Primary)	O15379
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500; FC-1:20; IP-1:20-1:50

Target information (O15379)

Gene Name	HDAC3
About Gene Symbol	HDAC3 (Histone Deacetylase 3) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC3 in 3T3, HeLa, rat Brain lysates using HDAC3 antibody.