

Hexokinase I Rabbit mAb


ABZ2985

Product Overview



Name	Hexokinase I Rabbit mAb
Accession(Primary)	P19367
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:20-1:50; FC-1:20

Target information (P19367)


Gene Name	HK1
About Gene Symbol	HK1 (Brain Form Hexokinase) is a protein kinase and signaling regulator involved in cellular metabolism. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Hexokinase I (green) in K562 using Hexokinase I antibody, and DAPI (blue).



Western blot analysis of Hexokinase I in K562, rat Brain, 3T3 lysates using Hexokinase I antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hexokinase I antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.