

Histone H2A Rabbit mAb

ABZ2987


Product Overview

Name	Histone H2A Rabbit mAb
Accession(Primary)	P04912
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Saccharomyces cerevisiae
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P04912)

Gene Name HTA2

About Gene Symbol



Western blot analysis of Histone H2A in K562, C6 lysates using Histone H2A antibody.

