

Histone H2A.X Rabbit mAb


ABZ2988

Product Overview


Name	Histone H2A.X Rabbit mAb
Accession(Primary)	P16104
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:200-1:1000; ICC/IF-1:200-1:1000; FC-1:20; IP-1:20-1:100

Target information (P16104)


Gene Name	H2AX
About Gene Symbol	




Immunocytochemistry analysis of Histone H2A.X (green) in HL-60 using Histone H2A.X antibody, and DAPI(blue).



Western blot analysis of Histone H2A.X in PC-12 lysates using Histone H2A.X antibody.



Western blot analysis of Histone H2A.X in 3T3, HeLa lysates using Histone H2A.X antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Histone H2A.X antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.