

Histone H2A.Z Rabbit mAb


ABZ2989

Product Overview

Name	Histone H2A.Z Rabbit mAb
Accession(Primary)	P0C0S5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (P0C0S5)

Gene Name	H2AZ1
About Gene Symbol	



Western blot analysis of Histone H2A.Z in K562, C6, 3T3 lysates using Histone H2A.Z antibody.