

HLA A Rabbit mAb


ABZ2991

Product Overview


Name	HLA A Rabbit mAb
Accession(Primary)	P04439
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:20

Target information (P04439)


Gene Name	HLA-A
About Gene Symbol	HLA-A (Major Histocompatibility Complex, Class I, A) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and



Western blot analysis of HLAA in K562, HeLa lysates using HLA A antibody.



Immunocytochemistry analysis of HLA A(green) in HeLa using HLA A antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HLA A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.