

HLA C Rabbit mAb


ABZ2992

Product Overview

Name	HLA C Rabbit mAb
Accession(Primary)	P10321
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:50-1:100

Target information (P10321)

Gene Name	HLA-C
About Gene Symbol	HLA-C (Major Histocompatibility Complex, Class I, C) is a membrane receptor involved in adaptive immunity. It is an important research target in autoimmune disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HLA-C in HeLa, 293T, Jurkat, 3T3, HepG2 lysates using HLA C antibody.