

hnRNP G T Rabbit mAb


ABZ2994

Product Overview

Name	hnRNP G T Rabbit mAb
Accession(Primary)	O75526
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; IF-1:10-1:100

Target information (O75526)

Gene Name	RBMXL2
About Gene Symbol	



Western blot analysis of hnRNP GT in K562, rat Brain, C6, 3T3, HeLa lysates using hnRNP G T antibody.