


hnRNP Q Rabbit mAb

ABZ2995


Product Overview

Name	hnRNP Q Rabbit mAb
Accession(Primary)	O60506
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (O60506)

Gene Name SYNCRIP**About Gene Symbol**


Western blot analysis of hnRNP Q in mouse intestinal, 3T3 lysates using hnRNP Q antibody.



Western blot analysis of hnRNP Q in K562 lysates using hnRNP Q antibody.



Immunocytochemistry analysis of hnRNP Q (green) in HeLa using hnRNP Q antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using hnRNP Q antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.