

HP1 gamma Rabbit mAb


ABZ2996

Product Overview


Name	HP1 gamma Rabbit mAb
Accession(Primary)	P23198
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:500-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (P23198)

Gene Name	Cbx3
About Gene Symbol	



Western blot analysis of HP1 gamma in K562, C6, 3T3 lysates using HP1 gamma antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HP1 gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.