


## HP1 gamma Rabbit mAb


ABZ2997

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | HP1 gamma Rabbit mAb                       |
| <b>Accession(Primary)</b>   | Q13185                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 21 kDa; Observed MW: 21 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:500-1:2000       |

### Target information (Q13185)

**Gene Name** CBX3**About Gene Symbol**



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HP1 gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.