

Hsc70 Rabbit mAb


ABZ2998

Product Overview


Name	Hsc70 Rabbit mAb
Accession(Primary)	P11142
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:10-1:50

Target information (P11142)


Gene Name	HSPA8
About Gene Symbol	HSPA8 (Heat Shock Protein Family A (Hsp70) Member 8) is a molecular chaperone involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mecha




Western blot analysis of Hsc70 in mouse heart, 3T3 lysates using Hsc70 antibody.



Western blot analysis of Hsc70 in HeLa lysates using Hsc70 antibody.



Western blot analysis of Hsc70 in Jurkat, rat Brain, C6, CHO-K1, Hela lysates using Hsc70 antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Hsc70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.