

HSD17B1 Rabbit mAb


ABZ2999

Product Overview

Name	HSD17B1 Rabbit mAb
Accession(Primary)	P14061
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:100-1:200; IP-1:20-1:50

Target information (P14061)

Gene Name	HSD17B1
About Gene Symbol	



Western blot analysis of HSD17B1 in Human placenta lysates using HSD17B1 antibody.