

Hsp60 Rabbit mAb


ABZ3001

Product Overview

Name	Hsp60 Rabbit mAb
Accession(Primary)	P10809
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (P10809)

Gene Name	HSPD1
About Gene Symbol	HSPD1 (Heat Shock Protein Family D (Hsp60) Member 1) is a molecular chaperone involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mecha



Western blot analysis of Hsp60 in rat Brain, 3T3
lysates using Hsp60 antibody.