

Hsp70 Rabbit mAb


ABZ3002

Product Overview


Name	Hsp70 Rabbit mAb
Accession(Primary)	P0DMV8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:10000; IHC-1:100-1:2000; ICC/IF-1:200-1:1000; FC-1:20

Target information (P0DMV8)


Gene Name	HSPA1A
About Gene Symbol	



Western blot analysis of Hsp70 in K562, rat Brain, C6, 3T3 lysates using Hsp70 antibody.



Immunocytochemistry analysis of Hsp70(green) in Hela using Hsp70 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.