

HtrA2 / Omi Rabbit mAb


ABZ3004

Product Overview

Name	HtrA2 / Omi Rabbit mAb
Accession(Primary)	O43464
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:50-1:100

Target information (O43464)

Gene Name	HTRA2
About Gene Symbol	



Western blot analysis of HtrA2/Omi in 3T3, HeLa lysates using HtrA2 antibody.