

ICAM1 Rabbit mAb


ABZ3005

Product Overview


Name	ICAM1 Rabbit mAb
Accession(Primary)	P05362
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 58 kDa; Observed MW: 89,92 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:10-1:100

Target information (P05362)

Gene Name	ICAM1
About Gene Symbol	



Western blot analysis of ICAM1 in K562, HeLa lysates using ICAM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ICAM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.