

ICOS Ligand/ICOSL Rabbit mAb

ABZ3006


Product Overview

Name	ICOS Ligand/ICOSL Rabbit mAb
Accession(Primary)	O75144
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (O75144)

Gene Name ICOSLG

About Gene Symbol



Western blot analysis of ICOS Ligand in K562 lysates using ICOS Ligand antibody.

