

IFITM1 Rabbit mAb


ABZ3007

Product Overview



Name	IFITM1 Rabbit mAb
Accession(Primary)	P13164
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50

Target information (P13164)


Gene Name	IFITM1
About Gene Symbol	



Immunocytochemistry analysis of IFITM1 (green) in K562 using IFITM1 antibody, and DAPI (blue).



Western blot analysis of IFITM1 in K562 lysates using IFITM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using IFITM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.