

IGFBP2 Rabbit mAb


ABZ3009

Product Overview

Name	IGFBP2 Rabbit mAb
Accession(Primary)	P18065
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (P18065)

Gene Name	IGFBP2
About Gene Symbol	



Western blot analysis of IGFBP2 in T47D lysates
using IGFBP2 antibody.