

IGFBP6 Rabbit mAb


ABZ3010

Product Overview

Name	IGFBP6 Rabbit mAb
Accession(Primary)	P24592
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:200-1:1000

Target information (P24592)

Gene Name	IGFBP6
About Gene Symbol	



Western blot analysis of IGFBP6 in K562, HeLa
lysates using IGFBP6 antibody.