

IL-1RA Protein Rabbit mAb


ABZ3013

Product Overview


Name	IL-1RA Protein Rabbit mAb
Accession(Primary)	P18510
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-1:50-1:100; ICC/IF-1:50-1:100

Target information (P18510)


Gene Name	IL1RN
About Gene Symbol	



Immunocytochemistry analysis of IL 1 Receptor Antagonist Protein (green) in HEPG2 using IL 1 Receptor Antagonist Protein antibody, and DAPI (blue).



Western blot analysis of IL-1RA in K562, 3T3, HeLa lysates using IL-1RA antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IL-1RA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.