

ILF3 Rabbit mAb


ABZ3014

Product Overview


Name	ILF3 Rabbit mAb
Accession(Primary)	Q12906
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (Q12906)


Gene Name	ILF3
About Gene Symbol	



Western blot analysis of ILF3/NF90 in Jurkat, C6, HeLa lysates using ILF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ILF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ILF3 (green) in
hela using ILF3 antibody, and DAPI(blue)