

IMP3 Rabbit mAb


ABZ3015

Product Overview


Name	IMP3 Rabbit mAb
Accession(Primary)	O00425
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20; IP-1:10-1:100

Target information (O00425)

Gene Name	IGF2BP3
About Gene Symbol	IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) is a cytokine and immune signaling protein involved in adaptive immunity and translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of IMP3 in K562, rat Brain, C6, 3T3, HeLa lysates using IMP3 antibody.



Immunocytochemistry analysis of IMP3(green) in HeLa using IMP3 antibody, and DAPI(blue)