

Indoleamine 2,3-dioxygenase Rabbit mAb


ABZ3016

Product Overview

Name	Indoleamine 2,3-dioxygenase Rabbit mAb
Accession(Primary)	P14902
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:500-1:2000; FC-1:200-1:500; IP-1:20-1:50

Target information (P14902)

Gene Name	IDO1
About Gene Symbol	



Western blot analysis of Indoleamine 2,
3-dioxygenase in Human placenta, SKOV-3 lysates
using Indoleamine 2, 3-dioxygenase antibody.