


Inhibin beta A Rabbit mAb

ABZ3017

Product Overview

Name	Inhibin beta A Rabbit mAb
Accession(Primary)	P08476
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P08476)

Gene Name INHBA**About Gene Symbol**

Western blot analysis of Inhibin beta A in rat Brain, C6, CHO-K1, HeLa lysates using Inhibin beta A antibody.

