

INPP4B Rabbit mAb


ABZ3018

Product Overview


Name	INPP4B Rabbit mAb
Accession(Primary)	O15327
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 105 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (O15327)

Gene Name	INPP4B
About Gene Symbol	



Western blot analysis of INPP4B in Hacat lysates using INPP4B antibody.



Western blot analysis of INPP4B in rat Brain, C6 lysates using INPP4B antibody.