

Integrin beta 3 Rabbit mAb


ABZ3019

Product Overview


Name	Integrin beta 3 Rabbit mAb
Accession(Primary)	P05106
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000

Target information (P05106)

Gene Name	ITGB3
About Gene Symbol	ITGB3 (Integrin Subunit Beta 3) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin beta 3 in U87-MG lysates using Integrin beta 3 antibody.



Western blot analysis of Integrin beta 3 in HeLa lysates using Integrin beta 3 antibody.