

Integrin beta 5 Rabbit mAb


ABZ3020

Product Overview

Name	Integrin beta 5 Rabbit mAb
Accession(Primary)	P18084
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P18084)

Gene Name	ITGB5
About Gene Symbol	ITGB5 (Integrin Subunit Beta 5) is an integrin cell adhesion receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin beta 5 in C6, CHO-K1, HeLa lysates using Integrin beta 5 antibody.