

ISG15 Rabbit mAb


ABZ3022

Product Overview

Name	ISG15 Rabbit mAb
Accession(Primary)	P05161
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (P05161)

Gene Name	ISG15
About Gene Symbol	ISG15 (ISG15 Ubiquitin Like Modifier) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling and ubiquitin-proteasome pathway. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ISG15 in C6, 3T3, HeLa
lysates using ISG15 antibody.