

Integrin beta 6 Rabbit mAb

ABZ3023


Product Overview

Name	Integrin beta 6 Rabbit mAb
Accession(Primary)	P18564
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P18564)

Gene Name ITGB6

About Gene Symbol



Western blot analysis of ITGB6 in rat Brain, C6, 3T3, HeLa lysates using Integrin beta 6 antibody.

