

KCTD3 Rabbit mAb


ABZ3029

Product Overview

Name	KCTD3 Rabbit mAb
Accession(Primary)	Q9Y597
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9Y597)

Gene Name	KCTD3
About Gene Symbol	KCTD3 (Potassium Channel Tetramerization Domain Containing 3) is a membrane transporter involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of KCTD3 in K562 lysates
using KCTD3 antibody.