

KGF/FGF-7 Rabbit mAb

ABZ3030


Product Overview

Name	KGF/FGF-7 Rabbit mAb
Accession(Primary)	P21781
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P21781)

Gene Name FGF7

About Gene Symbol



Western blot analysis of KGF in rat Brain, C6, 3T3, HeLa lysates using KGF antibody.

