


KLLN Rabbit mAb

ABZ3031

Product Overview

Name	KLLN Rabbit mAb
Accession(Primary)	B2CW77
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (B2CW77)

Gene Name KLLN**About Gene Symbol**

Western blot analysis of KLLN in K562, C6, 3T3, HeLa lysates using Killin antibody.

