

Lactate Dehydrogenase C Rabbit mAb


ABZ3033

Product Overview

Name	Lactate Dehydrogenase C Rabbit mAb
Accession(Primary)	P07864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P07864)

Gene Name	LDHC
About Gene Symbol	LDHC (Lactate Dehydrogenase C) is a metabolic dehydrogenase involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Lactate Dehydrogenase C in HeLa lysates using Lactate Dehydrogenase C antibody.