

Laminin gamma 1 Rabbit mAb


ABZ3034

Product Overview


Name	Laminin gamma 1 Rabbit mAb
Accession(Primary)	P11047
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 178 kDa; Observed MW: 220-250 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; FC-1:20-1:50; IP-1:20-1:50

Target information (P11047)

Gene Name	LAMC1
About Gene Symbol	LAMC1 (Laminin Subunit Gamma 1) is an extracellular matrix laminin protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Laminin gamma 1 in C6, mouse brain lysates using Laminin gamma 1 antibody.



Western blot analysis of Laminin gamma 1 in HeLa, MCF-7, A549, HCT116 lysates using Laminin gamma 1 antibody.