

LAMP1 Rabbit mAb


ABZ3035

Product Overview


Name	LAMP1 Rabbit mAb
Accession(Primary)	P11279
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 90-120 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500

Target information (P11279)

Gene Name	LAMP1
About Gene Symbol	LAMP1 (Lysosomal Associated Membrane Protein 1) is a cell surface marker involved in adaptive immunity and vesicular trafficking. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of LAMP1 in HeLa lysates using LAMP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LAMP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.