

LAMP2 Rabbit mAb


ABZ3036

Product Overview


Name	LAMP2 Rabbit mAb
Accession(Primary)	P13473
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:20-1:100; FC-1:100-1:200; IP-1:20

Target information (P13473)

Gene Name	LAMP2
About Gene Symbol	LAMP2 (Lysosomal Associated Membrane Protein 2) is a cell surface marker involved in adaptive immunity and vesicular trafficking. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of LAMP2 in 3T3 lysates using LAMP2A antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using LAMP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.