

Lck Rabbit mAb


ABZ3037

Product Overview


Name	Lck Rabbit mAb
Accession(Primary)	P06239
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P06239)


Gene Name	LCK
About Gene Symbol	




Immunocytochemistry analysis of Lck (green) in Jurkat using Lck antibody, and DAPI (blue).



Western blot analysis of Lck in rat spleen lysates using Lck antibody.



Western blot analysis of Lck in Jurkat lysates
using Lck antibody



Immunohistochemistry analysis of paraffin-
embedded Human tonsil using Lck antibody. High-
pressure and temperature Sodium Citrate pH 6.0
was used for antigen retrieval.