

Leukotriene A4 Hydrolase Rabbit mAb


ABZ3038

Product Overview


Name	Leukotriene A4 Hydrolase Rabbit mAb
Accession(Primary)	P09960
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:100-1:1000

Target information (P09960)

Gene Name	LTA4H
About Gene Symbol	LTA4H (Leukotriene A4 Hydrolase) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Leukotriene A4 hydrolase/LTA4H in K562, C6, 3T3, HeLa lysates using Leukotriene A4 Hydrolase antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Leukotriene A4 hydrolase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.