

LGALS14 Rabbit mAb


ABZ3039

Product Overview



Name	LGALS14 Rabbit mAb
Accession(Primary)	Q8TCE9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (Q8TCE9)

Gene Name	LGALS14
About Gene Symbol	



Immunocytochemistry analysis of LGALS14 (green) in K562 using LGALS14 antibody, and DAPI (blue).



Western blot analysis of LGALS14 in C6 lysates using LGALS14 antibody.