

liver FABP Rabbit mAb


ABZ3040

Product Overview

Name	liver FABP Rabbit mAb
Accession(Primary)	P07148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100

Target information (P07148)

Gene Name	FABP1
About Gene Symbol	



Western blot analysis of Liver FABP in mouse liver lysates using FABP1 antibody.