

LRG1 Rabbit mAb


ABZ3041

Product Overview

Name	LRG1 Rabbit mAb
Accession(Primary)	P02750
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P02750)

Gene Name	LRG1
About Gene Symbol	LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of LRG1 in Human serum
lysates using LRG1 antibody.