

mH2A1 Rabbit mAb


ABZ3042

Product Overview


Name	mH2A1 Rabbit mAb
Accession(Primary)	O75367
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500

Target information (O75367)

Gene Name	MACROH2A1
About Gene Symbol	MACROH2A1 (MacroH2A.1 Histone) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and adaptive immunity. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of macroH2A.1 in rat Brain, C6, 3T3 lysates using Macro H2A.1 antibody.



Immunocytochemistry analysis of macroH2A.1(green) in Hela using macroH2A.1 antibody, and DAPI(blue)