

MAGE3 Rabbit mAb

ABZ3043


Product Overview

Name	MAGE3 Rabbit mAb
Accession(Primary)	P43357
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P43357)

Gene Name MAGEA3

About Gene Symbol



Western blot analysis of MAGEA3 in HeLa lysates using MAGEA3 antibody.

